








EDGEFILE® X7.06

Specifications Overview

	Tip Size	Taper	Max Diameter	16mm	13mm	9mm	6mm	3mm	1mm	Torque g/cm	Speed RPM
	17	0.06	1.0	1.0	0.95	0.71	0.53	0.35	0.23	250 to 410	300 to 500
	20	0.06	1.0	1.0	0.98	0.74	0.56	0.38	0.26	250 to 410	300 to 500
	25	0.06	1.0	1.0	1.0	0.79	0.61	0.43	0.31	250 to 410	300 to 500
	30	0.06	1.0	1.0	1.0	0.84	0.66	0.48	0.36	250 to 410	300 to 500
	35	0.06	1.0	1.0	1.0	0.89	0.71	0.53	0.41	250 to 410	300 to 500
	40	0.06	1.0	1.0	1.0	0.94	0.76	0.58	0.46	250 to 410	300 to 500
	45	0.06	1.0	1.0	1.0	0.99	0.81	0.63	0.51	250 to 410	300 to 500