





Specifications Overview

Reverse reciprocating 30° clockwise and 150° counterclockwise	Tip Size	Max Diameter	16mm	13mm	9mm	6mm	3mm	1mm	Torque g/cm N/cm
Small 	0.20	1.20	1.20	0.87	0.72	0.54	0.37	0.25	408 – 530 g/cm 4,0 – 5,2 N/cm
Primary 	0.25	1.20	1.20	0.98	0.82	0.61	0.43	0.30	408 – 530 g/cm 4,0 – 5,2 N/cm
Medium 	0.35	1.20	1.20	0.95	0.82	0.63	0.49	0.40	408 – 530 g/cm 4,0 – 5,2 N/cm
Large 	0.45	1.20	1.20	0.92	0.82	0.66	0.56	0.50	408 – 530 g/cm 4,0 – 5,2 N/cm