








# EDGE X7 UTOPIA™

## SPECIFICATIONS OVERVIEW .06

	Tip Size	Taper	Max Diameter	16mm	13mm	9mm	6mm	3mm	1mm	Torque [N/cm]	Speed [RPM]
	17	0.06	1.0	1.0	0.95	0.71	0.53	0.35	0.23	4.1	500
	20	0.06	1.0	1.0	0.98	0.74	0.56	0.38	0.26	4.1	500
	25	0.06	1.0	1.0	1.0	0.79	0.61	0.43	0.31	4.1	500
	30	0.06	1.0	1.0	1.0	0.84	0.66	0.48	0.36	4.1	500
	35	0.06	1.0	1.0	1.0	0.89	0.71	0.53	0.41	4.1	500
	40	0.06	1.0	1.0	1.0	0.94	0.76	0.58	0.46	4.1	500
	45	0.06	1.0	1.0	1.0	0.99	0.81	0.63	0.51	4.1	500