







# EDGE TAPER BLAZE UTOPIA™

## SPECIFICATIONS OVERVIEW

	Tip Size	Taper	Max Diameter	16mm	13mm	9mm	6mm	3mm	1mm	Torque [g/cm]	Speed [RPM]
 SX	19	04v	1.2	1.2	1.09	1.09	0.58	0.32	0.22	520	300
 S1	18	03v	1.2	1.2	1.06	0.65	0.41	0.26	0.19	520	300
 S2	20	04v	1.2	1.2	1.05	0.69	0.51	0.33	0.24	153	300
 F1	20	07v	1.2	1.12	0.96	0.74	0.57	0.41	0.27	153	300
 F2	25	08v	1.2	1.09	0.97	0.81	0.66	0.49	0.33	316	300
 F3	30	09v	1.2	1.03	1.03	0.87	0.75	0.57	0.39	316	300

©EdgeEndo 2023/04