





# EDGE ONE R UTOPIA™

## SPECIFICATIONS OVERVIEW

	Tip Size	Taper	Max Diameter	16mm	13mm	9mm	6mm	3mm	1mm	Torque [4N/cm]	Speed [RPM]
	20	08v	1.2	1.0	0.88	0.72	0.60	0.44	0.28	170°/60°	400
	25	08v	1.2	1.5	0.93	0.77	0.65	0.49	0.33	170°/60°	400
	40	06v	1.2	1.1	0.98	0.82	0.70	0.58	0.46	170°/60°	400
	50	05v	1.2	1.17	1.05	0.89	0.77	0.65	0.55	170°/60°	400

©EdgeEndo 2023/04