



# Treating with precision.

## ZEISS EyeMag Medical Loupes


**Customer Presentation**

ZEISS Medical Technology



# ZEISS EyeMag Medical Loupes

## Two Centuries of Precision



### Visionary breakthrough

The visionary development is followed by the patent application on September 10, 1909.



### Modern vision with Galilean design

Development of Galilean-style loupes in a small design with 2.0x magnification.



### Monocular masterpiece

Introduction of the Keplerian monocular and binocular telescope with 6.0x magnification.



### Shaping a new century

Relocation of medical loupes from ZEISS Vision to Medical Technology, leading to the introduction of the Smart, Pro F, and Pro S models.




### Next level of illumination

Sleeker, lighter, and more portable than previous ZEISS models – evolution to modernity with EyeMag Light II.



### What's next?

Curious about the next evolution of the EyeMag® Medical Loupes from ZEISS?



### **Advanced perspective**

With the addition of the EyeMag Classic and EyeMag Smart TTL, the portfolio now includes individualized and through-the-lens (TTL) loupes for enhanced performance.

# ZEISS EyeMag Medical Loupes


As individual as your needs



**Redesigned Optics**



**Personal Choice**



**Intuitive Handling**

# ZEISS EyeMag Medical Loupes

## Redesigned Optics



Seeing tiny anatomical structures in detail is essential for a successful treatment. Clear and sharp images are what the wearers expect.

The high-quality ZEISS optics enable precise visual presentation of minute details. ZEISS EyeMag medical loupes provide excellent optical performance and a wide range of magnification to cover all your needs.

# ZEISS EyeMag Medical Loupes

## Personal Choice



Choose from a wide variety of configurations. The large field of view and depth of field of the Galilean optics enable easy visualization of the entire treatment area. The Kepler optics provide high-contrast stereoscopic images for detail recognition and differentiation. Both optics are available as individualized or flip-up versions and can be paired with lightweight frames in a variety of different colors. You can also complement your medical loupe with the powerful lighting of the ZEISS EyeMag Light II.

# ZEISS EyeMag Medical Loupes

## Intuitive Handling



ZEISS medical loupes allow you to adapt to the respective situation. Simply put on your pair of individualized medical loupes and you are good to go. When using the flip-up medical loupes, the eyepieces can be repositioned precisely with a single adjustment or even flipped up entirely for eye contact with your patient when magnification is not needed.

# ZEISS EyeMag Medical Loupes

## Ergonomic benefits



01  
Ergonomic  
position



02  
Non-ergonomic  
position



### Did you know...

...that back and neck pain due to awkward positioning and posture is one of the most common challenges dentists are facing?

Due to their compact, portable, and ergonomic design, ZEISS EyeMag medical loupes are supporting an ergonomic working position, that provides maximum comfort in daily use for a wide range of procedures.

# ZEISS EyeMag Medical Loupes

The entire product portfolio



EyeMag Classic



EyeMag Smart TTL



EyeMag Light II



EyeMag Smart



EyeMag Pro F / S

# ZEISS EyeMag Medical Loupes

New frames and colors for individualized loupes



Nova



Polaris



Sirius



Felis



Hoshi



Realta



Find the perfect medical loupe using the [product configurator](#).

# ZEISS EyeMag Medical Loupes

Frames for flip-up and individualized loupes



STMS



Strator\*



Find the perfect medical loupe using the [product configurator](#).

\* Only available for the EyeMag Smart

# ZEISS EyeMag Medical Loupes

## Personalize your medical loupe



So many options, find the right one for you with the **product configurator!**

Your medical loupe is as personal as your pair of glasses – it has to provide great optical performance, while also suiting your personal style. To help you get an overview of the ZEISS EyeMag medical loupes product portfolio, we created the product configurator. Select the type of loupe you want to learn more about and configure the frame and color according to your preferences.

[Configure your medical loupe](#)

# ZEISS EyeMag Medical Loupes

## EyeMag Classic



- Individualized medical loupe, the optical systems are mounted on the carrier lenses surface
- With 2.0x or 2.5x magnification
- Available with titanium eyeglass frames in various colors
- Individual prescription values within system
- Attractively-priced entry model, ideal for users with basic magnification requirements or students

# ZEISS EyeMag Medical Loupes

## EyeMag Smart



- Flip-up
- With 2.5x magnification
- Available with titanium eyeglass frame or sports frame
- Ideal for users with basic magnification requirements and need for flexible angle change during usage

# ZEISS EyeMag Medical Loupes

## EyeMag Smart TTL



- Through the lens, the optical systems are mounted into the carrier lenses
- With 2.5x or 3.0x magnification
- Available with titanium eyeglass frames in various colors
- Individual prescription values within system
- Ideal for users with basic magnification requirements and who need large field of view

# ZEISS EyeMag Medical Loupes

## EyeMag Pro F



- Flip-up
- With magnification from 3.2x to 5.0x
- Available with titanium eyeglass frame
- Ideal for users with high magnification needs and require flexible angle change during usage

# ZEISS EyeMag Medical Loupes

## EyeMag Pro S



- Flip-up
- With magnification from 3.2x to 5.0x
- Available with a comfortable elastic strap
- Ideal for users with high magnification needs and require flexible angle change during usage

# ZEISS EyeMag Medical Loupes

## EyeMag Light II



- Combined with high-quality ZEISS optics, the powerful ZEISS EyeMag Light II provides crisp detail recognition and an unaltered view of tiny structures
- The almost coaxial light results in minimal shadow formation
- Two powerful rechargeable battery packs and a charging cradle allows for practically interruption-free work – 24 / 7
- Perfect fit for all ZEISS EyeMag medical loupes



Seeing beyond